

LABORATORY HOT PLATE

- Rugged, steel construction with solid cast iron plate on top.
- Heaters are designed and distributed for uniform surface temperature.
- With capillary thermostat for accurate temperature control.
- Electrically operated on 230 VAC, Single Phase.



ROUND MODEL

Model Code	Size (dia)		Rating
HHPC - 8	20 cm.	8 inch	1 kW
HHPC - 12	30 cm.	12 inch	1.5 kW



RECTANGULAR MODEL

Model Code	Size (L x B)		Rating
HHPC - 1012	25 x 30cm.	10 x 12 inch	1 kW
HHPC - 1218	30 x 45cm.	12 x 18 inch	1.5 kW

HARIDAS INSTRUMENTS & EQUIPMENTS CO.

(An ISO 9001:2008 Company)

Sector-23, Plot No.76, CIDCO Service Industrial Area, Turbhe, Navi Mumbai-400 705, Maharashtra, India.
 Tel : 91-022-27833046, 27835042 Fax : 91-022-27835042 email : kph@vsnl.com url : [http:// www.haridasinstruments.com](http://www.haridasinstruments.com)