


LABORATORY MUFFLE FURNACE



- Heavy-duty steel exterior construction.
- Interior of ceramic muffle.
- Thick ceramic blanket insulation provides excellent protection from heat loss.
- Temperature range upto 950°C.
- Choice of digital On/Off or PID temperature control.
- Electrically operated on 230 VAC, Single Phase.

Model Code	Muffle Size (H x W x D)	Rating
HMF - 4	10x10x23cm	1.6 kW
HMF - 5	13x13x25cm	2.6 kW
HMF - 6	15x15x30cm	3.3 kW

HARIDAS INSTRUMENTS & EQUIPMENTS CO.
(An ISO 9001:2008 Company)

Sector-23, Plot No.76, CIDCO Service Industrial Area, Turbhe, Navi Mumbai-400 705, Maharashtra, India.
 Tel : 91-022-27833046, 27835042 Fax : 91-022-27835042 email : kph@vsnl.com url : http:// www.haridasinstruments.com