

LABORATORY WATER BATH

- Double walled with interior & lid of stainless steel.
- Glasswool insulation ensures minimal heat loss.
- Provided with concentric rings to accommodate various sizes of glass flasks.
- Choice of energy regulator or thermostatic temperature control.
- Electrically operated on 230 VAC, Single Phase.



ROUND MODEL

Model Code	Flask Size	Rating
HWRF - 2	2 Lt.	1 kW
HWRF - 5	5 Lt.	1.5 kW



RECTANGULAR MODEL (THERMOSTATIC)

Model Code	Size (L x B x H)	No. of Holes	Rating
HWB - 6	35 x 25 x 10cm	6	1 kW
HWB - 12	40 x 30 x 10cm	12	2 kW

HARIDAS INSTRUMENTS & EQUIPMENTS CO.

(An ISO 9001:2008 Company)

Sector-23, Plot No.76, CIDCO Service Industrial Area, Turbhe, Navi Mumbai-400 705, Maharashtra, India.
 Tel : 91-022-27833046, 27835042 Fax : 91-022-27835042 email : kph@vsnl.com url : [http:// www.haridasinstruments.com](http://www.haridasinstruments.com)